

## CLINICAL 360°

### TECH BRIEF

#### APPLYING

## MIDWEST STYLUS ATC *in your practice*

Compiled by **RICHARD PALMER**  
Information provided by **DENTSPLY**  
Midwest.

The lightweight, ergonomic **Midwest® Stylus™ ATC** highspeed air-driven handpiece delivers an exceptionally smooth, powerful and efficient cutting performance.

**THE NEED FOR SPEED:** The Speed-Sensing Intelligence (SSI) constantly maintains bur speed according to load and automatically optimizes power to virtually eliminate stalling and the need to make frequent adjustments.

**BELIEF OF SUSPENSION:** Superior Turbine Suspension (SPS) allows the handpiece to operate at speeds of 330,000 rpm without noticeable bur deflection or chatter for optimal control, precision and transference of power to the bur.

**GRACE UNDER PRESSURE:** ATC Technology adjusts air pressure to the handpiece depending on bur load to deliver optimal cutting performance when needed.

**ADDED PLUS:** The user-friendly design features a swivel with exceptional rotation capabilities, lightweight ergonomics, exceptional visibility, enhanced access to posterior surfaces, and low pitch and tone to minimize stress for clinician and patient. ●



ARE YOU READY  
TO APPLY THIS  
TECHNOLOGY IN  
YOUR PRACTICE?

**DENTSPLY Midwest**  
[StylusATC.com](http://StylusATC.com)  
800-989-8825